

Single Use Retrieval Basket V



Single Use Retrieval Basket V



# Revolutionary V-System Device Exchange is Now Combined with New Basket Designs, a Rotation Mechanism, and Wire Guided Access



## The Innovative V-System Design Lets You Proceed with Confidence and Efficiency


The V-System is a complete system that integrates Olympus endoscopes and EndoTherapy devices. The revolutionary V-System design offers the option of guidewire manipulation by the physician or the assistant, allows easier exchange of catheters, and enhances cannulation capability.



### V-Marking

Indicates when to raise and lower the V-Groove forceps elevator.

The exclusive V-Marking is located on the proximal side of the sheath. When this marking reaches the channel port on the scope's control section, it indicates that the device tip has reached the distal end of the scope and the V-Groove forceps elevator may be lowered. When withdrawing the device from the scope, the same marking indicates when to raise the elevator to lock the guidewire.



Unique Distal Eight-wire Basket is Ideal for Retrieval of Small Stones and with Four Proximal Wires to Facilitate Release of Captured Stones



Standard Four-wire Basket Meets Most Stone Retrieval Requirements in the Bile and Pancreatic Ducts


#### Wire guided ability for superior insertion capability

The FG-V411Q and FG-V412Q are designed to be wire guided with the wire passing only through the distal end of the basket. This makes it possible to insert the basket into a small incision on the papilla, and allows selective insertion into the intrahepatic bile duct. In addition, these baskets have the same expansion capability as our conventional models (FG-401Q/FG-22Q-1 etc) because the guidewire is only attached to the basket's distal tip.




#### Eight-wire basket designed for the capture and removal of small stones.

The FlowerBasketV has eight wires on the distal side and four wires on the proximal side. The eight-wire section allows the precise retrieval of small stones while the four-wire section makes it easier to release them.



#### Four-wire basket design for versatile retrieval

For versatile retrieval of a wide range of stones, the TetraCatchV features a standard four-wire retrieval basket configuration.



#### Rotation mechanism and bullet-shaped tip

The rotatable basket (FG-V401QR/FG-V402QR) allows the retrieval of small stones, especially floating stones or stones located in the lower bile duct above the papilla. The bullet-shaped tip ensures smooth insertion into the bile duct.



- Injection port for contrast media facilitates fluoroscopic visualization.
- Sterile single-use design for convenience both before and after procedures.

## Specifications

Model	Basket type	Rotatable/Guidewire type	Working length	Opening width	Minimum channel size	Compatible guidewire
FG-V401QR	8wire type	Rotatable type	1950mm	20mm	2.8mm	-
FG-V402QR	4wire type	Rotatable type	1950mm	22mm	2.8mm	-
FG-V411Q	8wire type	Guidewire type	1950mm	20mm	3.7mm	0.89mm(0.035inch)
FG-V412Q	4wire type	Guidewire type	1950mm	22mm	3.7mm	0.89mm(0.035inch)



### BML-110A-1

An emergency device, this is only to be used in the unlikely event of basket impaction due to an extremely hard calculus.

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

**OLYMPUS**<sup>®</sup>

www.olympus.com

**OLYMPUS CORPORATION**  
Shinjuku Monolith, 3-1 Nishi-Shinjuku 2-chome, Shinjuku-ku, Tokyo 163-0914, Japan  
**OLYMPUS EUROPA GMBH**  
Postfach 10 49 08, 20034 Hamburg / Wendenstrasse 14-18, 20097 Hamburg, Germany  
**OLYMPUS AMERICA INC.**  
2 Corporate Center Drive, Melville, N.Y. 11747-3157, U.S.A.  
**OLYMPUS LATIN AMERICA, INC.**  
6100 Blue Lagoon Drive, Suite 390 Miami, Florida 33126-2087, U.S.A  
**KEYMED LTD.**  
Keymed House, Stock Road, Southend-on-Sea, Essex SS2 5QH, England

**OLYMPUS SINGAPORE PTE LTD.**  
491B River Valley Road #12-01/04, Valley Point Office Tower, Singapore 248373  
**OLYMPUS HONG KONG AND CHINA LIMITED**  
Room 1520-1527, Ocean Centre, 5 Canton Road, Tsimshatsui, Kowloon, Hong Kong  
**OLYMPUS BEIJING REPRESENTATIVE OFFICE**  
Room 1401, East Tower, Gong Yuan No8 Royal Palace, No8 Gong Yuan Xie, JianGuoMenNei, DongCheng District, Beijing, 100005, China  
**OLYMPUS MOSCOW LIMITED LIABILITY COMPANY**  
117071, Moscow, Malaya Kaluzhskaya 19, bld. 1, fl.2, Russia  
**OLYMPUS AUSTRALIA PTY LTD.**  
31 Gilby Road, Mount Waverley, VIC., 3149, Australia

Printed in Japan F1091SB-3-0804