



Your Vision, Our Future

February 24, 2011

No Microsoft Patch test executes for Windows 2000 starting month of March

The following list of Microsoft **Windows 2000** updates have been tested and approved for **EndoWorks 7** compatibility. Updates marked in red have associated Microsoft Security Bulletins. Prior to applying Windows Updates, make sure your system is current with EndoWorks patches. For questions regarding applying Microsoft Updates or EndoWorks Patches, please contact Olympus service.

EndoWorks does not support Internet Explorer 7.

Windows Updates	Description	Tested	Pass/Fail	Date
KB890830	Windows Malicious Software Removal Tool - June 2010	Yes	Pass	6/25/2010
KB975562	Security Update for DirectX 9 for Windows 2000	Yes	Pass	6/25/2010
KB979482	Security Update for Windows 2000	Yes	Pass	6/25/2010
KB979559	Security Update for Windows 2000	Yes	Pass	6/25/2010
KB979906	Microsoft .NET Framework 1.1 SP1 Security Update for Windows 2000	Yes	Pass	6/25/2010
KB980195	Cumulative Security Update for ActiveX Killbits for Windows 2000	Yes	Pass	6/25/2010
KB980218	Security Update for Windows 2000	Yes	Pass	6/25/2010
KB982381	Cumulative Security Update for Internet Explorer 6 Service Pack 1	Yes	Pass	6/25/2010
KB890830	Windows Malicious Software Removal Tool - May 2010	Yes	Pass	5/26/2010
KB978542	Security Update for outlook Express 6.0 for Windows 2000	Yes	Pass	5/26/2010
KB890830	Windows Malicious Software Removal Tool - April 2010	Yes	Pass	4/26/2010
KB977816	Security Update for Windows 2000	Yes	Pass	4/26/2010
KB978601	Security Update for Windows 2000	Yes	Pass	4/26/2010
KB979309	Security Update for Windows 2000	Yes	Pass	4/26/2010
KB979402	Security Update for Windows Media Player 9 for Windows 2000	Yes	Pass	4/26/2010
KB979683	Security Update for Windows 2000	Yes	Pass	4/26/2010
KB980182	Cumulative Security Update for Internet Explorer 6 Service Pack 1	Yes	Pass	4/26/2010
KB980232	Security Update for Windows 2000	Yes	Pass	4/26/2010
KB981350	Security Update for Windows 2000	Yes	Pass	4/26/2010
KB890830	Windows Malicious Software Removal Tool - March 2010	Yes	Pass	3/16/2010
KB890830	Windows Malicious Software Removal Tool - February 2010	Yes	Pass	2/22/2010
KB971468	Security Update for Windows 2000	Yes	Pass	2/22/2010
KB975560	Security Update for DirectX 9 for Windows 2000	Yes	Pass	2/22/2010
KB977165	Security Update for Windows 2000	Yes	Pass	2/22/2010
KB977914	Security Update for Windows 2000	Yes	Pass	2/22/2010
KB978037	Security Update for Windows 2000	Yes	Pass	2/22/2010
KB978251	Security Update for Windows 2000	Yes	Pass	2/22/2010
KB978262	Cumulative Security Update for ActiveX Killbits for Windows 2000	Yes	Pass	2/22/2010
KB978706	Security Update for Windows 2000	Yes	Pass	2/22/2010
KB890830	Windows Malicious Software Removal Tool - December 2009	Yes	Pass	12/22/2009
KB911565	Security Update for Windows Media Player 9	Yes	Pass	12/22/2009
KB923694	Cumulative Security Update for Outlook Express 6 Service Pack 1	Yes	Pass	12/22/2009
KB924883	Security Update for Word Viewer 2003	Yes	Pass	12/22/2009
KB951748	Security Update for Windows 2000	Yes	Pass	12/22/2009
KB955759	Update for Windows 2000	Yes	Pass	12/22/2009
KB973904	Security Update for Windows 2000	Yes	Pass	12/22/2009
KB974318	Security Update for Windows 2000	Yes	Pass	12/22/2009
KB974392	Security Update for Windows 2000	Yes	Pass	12/22/2009
KB976138	Update for DirectX 9 for Windows 2000	Yes	Pass	12/22/2009
KB976325	Cumulative Security Update for Internet Explorer 6 Service Pack 1	Yes	Pass	12/22/2009
KB969947	Security Update for Windows 2000	Yes	Pass	12/10/2009
KB890830	Windows Malicious Software Removal Tool - September 2009	Yes	Pass	10/24/2009
KB953155	Security Update for Windows 2000	Yes	Pass	10/24/2009
KB958869	Security Update for Internet Explorer 6 for Windows 2000	Yes	Pass	10/24/2009
KB969059	Security Update for Windows 2000	Yes	Pass	10/24/2009
KB971486	Security Update for Windows 2000	Yes	Pass	10/24/2009

KB973525	Cumulative Security Update for ActiveX Killbits for Windows 2000	Yes	Pass	10/24/2009
KB974112	Security Update for Windows 2000	Yes	Pass	10/24/2009
KB974455	Cumulative Security Update for Internet Explorer 6 Service Pack 1	Yes	Pass	10/24/2009
KB974571	Security Update for Windows 2000	Yes	Pass	10/24/2009
KB890830	Windows Malicious Software Removal Tool - September 2009	Yes	Pass	9/10/2009
KB956844	Security Update for Windows 2000	Yes	Pass	9/10/2009
KB968816	Security Update for Windows Media Format Runtime 9 for windows 2000	Yes	Pass	9/10/2009
KB971961	Security Update for Windows 2000	Yes	Pass	9/10/2009
KB890830	Windows Malicious Software Removal Tool - August 2009	Yes	Pass	8/12/2009
KB958470	Security Update for Windows 2000	Yes	Pass	8/12/2009
KB960859	Security Update for Windows 2000	Yes	Pass	8/12/2009
KB961371	Security Update for Windows 2000	Yes	Pass	8/12/2009
KB971557	Security Update for Windows 2000	Yes	Pass	8/12/2009
KB972260	Cumulative Security Update for Internet Explorer 6 Service Pack 1	Yes	Pass	8/12/2009
KB973354	Security Update for outlook Express 6.0 for Windows 2000	Yes	Pass	8/12/2009
KB973507	Security Update for Windows 2000	Yes	Pass	8/12/2009
KB890830	Windows Malicious Software Removal Tool - July 2009	Yes	Pass	7/15/2009
KB961371	Security Update for Windows 2000	Yes	Pass	7/15/2009
KB971633	Security Update for DirectX 9 for Windows 2000	Yes	Pass	7/15/2009
KB973346	Cumulative Security Update for ActiveX Killbits for Windows 2000	Yes	Pass	7/15/2009
KB890830	Windows Malicious Software Removal Tool - June 2009	Yes	Pass	6/25/2009
KB961501	Security Update for Windows 2000	Yes	Pass	6/25/2009
KB968537	Security Update for Windows 2000	Yes	Pass	6/25/2009
KB969897	Cumulative Security Update for Internet Explorer 5.01 Service Pack 4	Yes	Pass	6/25/2009
KB969897	Cumulative Security Update for Internet Explorer 6 Service Pack 1	Yes	Pass	6/25/2009
KB969898	Update Rollup for ActiveX Killbits for Windows 2000	Yes	Pass	6/25/2009
KB970238	Security Update for Windows 2000	Yes	Pass	6/25/2009
KB890830	Windows Malicious Software Removal Tool - May 2009	Yes	Pass	5/16/2009
KB890830	Windows Malicious Software Removal Tool - April 2009	Yes	Pass	4/14/2009
KB923561	Security Update for Windows 2000	Yes	Pass	4/14/2009
KB952004	Security Update for Windows 2000	Yes	Pass	4/14/2009
KB959426	Security Update for Windows 2000	Yes	Pass	4/14/2009
KB960803	Security Update for Windows 2000	Yes	Pass	4/14/2009
KB961373	Security Update for DirectX 9 for Windows 2000	Yes	Pass	4/14/2009
KB963027	Cumulative Security Update for Internet Explorer 6 Service Pack 1	Yes	Pass	4/14/2009
KB890830	Windows Malicious Software Removal Tool - March 2009	Yes	Pass	3/24/2009
KB958690	Security Update for Windows 2000	Yes	Pass	3/24/2009
KB960225	Security Update for Windows 2000	Yes	Pass	3/24/2009
KB967715	Update for Windows 2000	Yes	Pass	3/24/2009
KB890830	Windows Malicious Software Removal Tool - February 2009	Yes	Pass	2/17/2009
KB960715	Update Rollup for ActiveX Killbits for Windows 2000	Yes	Pass	2/17/2009
KB890830	Windows Malicious Software Removal Tool - January 2009	Yes	Pass	1/14/2009
KB960714	Security Update for Internet Explorer 6 Service Pack 1	Yes	Pass	12/17/2008
KB890830	Windows Malicious Software Removal Tool - December 2008	Yes	Pass	12/12/2008
KB952069	Security Update for Windows 2000	Yes	Pass	12/12/2008
KB954600	Security Update for Windows 2000	Yes	Pass	12/12/2008
KB956802	Security Update for Windows 2000	Yes	Pass	12/12/2008
KB958215	Cumulative Security Update for Internet Explorer 6 Service Pack 1	Yes	Pass	12/12/2008
KB890830	Windows Malicious Software Removal Tool - November 2008	Yes	Pass	11/13/2008
KB955069	Security Update for Windows 2000	Yes	Pass	11/13/2008
KB957097	Security Update for Windows 2000	Yes	Pass	11/13/2008
KB958644	Security Update for Windows 2000 (Critical Update)	Yes	Pass	10/27/2008
KB890830	Windows Malicious Software Removal Tool - October 2008	Yes	Pass	10/16/2008
KB954211	Security Update for Windows 2000	Yes	Pass	10/16/2008
KB956390	Cumulative Security Update for Internet Explorer 6 Service Pack 1	Yes	Pass	10/16/2008
KB956391	Cumulative Security Update for ActiveX Killbits for Windows 2000	Yes	Pass	10/16/2008
KB957095	Security Update for Windows 2000	Yes	Pass	10/16/2008
KB890830	Windows Malicious Software Removal Tool - September 2008	Yes	Pass	9/29/2008
KB938464	Security Update for Internet Explorer 6 for Windows 2000	Yes	Pass	9/29/2008
KB947742	Security Update for Microsoft .NET Framework, Version 1.1, Service Pack 1	Yes	Pass	9/29/2008
KB890830	Windows Malicious Software Removal Tool - August 2008	Yes	Pass	8/21/2008
KB950974	Security Update for Windows 2000	Yes	Pass	8/21/2008
KB952954	Security Update for Windows 2000	Yes	Pass	8/21/2008
KB953838	Cumulative Security Update for Internet Explorer 6 Service Pack 1	Yes	Pass	8/21/2008

KB953839	Cumulative Security Update for ActiveX Killbits for Windows 2000	Yes	Pass	8/21/2008
KB890830	Windows Malicious Software Removal Tool - July 2008	Yes	Pass	8/11/2008
KB951748	Security Update for Windows 2000	Yes	Pass	8/11/2008
KB890830	Windows Malicious Software Removal Tool - June 2008	Yes	Pass	6/18/2008
KB950759	Cumulative Security Update for Internet Explorer 6 Service Pack 1	Yes	Pass	6/18/2008
KB950760	Cumulative Security Update for ActiveX Killbits for Windows 2000	Yes	Pass	6/18/2008
KB951698	Security Update for DirectX 9 for Windows 2000	Yes	Pass	6/18/2008
KB890830	Windows Malicious Software Removal Tool - May 2008	Yes	Pass	5/30/2008
KB950749	Security Update for Windows 2000	Yes	Pass	5/30/2008
KB890830	Windows Malicious Software Removal Tool - April 2008	Yes	Pass	4/22/2008
KB941693	Security Update for Windows 2000	Yes	Pass	4/22/2008
KB944338	Security Update for Windows 2000	Yes	Pass	4/22/2008
KB945553	Security Update for Windows 2000	Yes	Pass	4/22/2008
KB947864	Cumulative Security Update for Internet Explorer 6 Service Pack 1	Yes	Pass	4/22/2008
KB948590	Security Update for Windows 2000	Yes	Pass	4/22/2008
KB948881	Security Update for ActiveX Killbits for Internet Explorer 6 for Windows 2000	Yes	Pass	4/22/2008
KB890830	Windows Malicious Software Removal Tool - March 2008	Yes	Pass	3/18/2008
KB890830	Windows Malicious Software Removal Tool - February 2008	Yes	Pass	2/22/2008
KB943055	Security Update for Windows 2000	Yes	Pass	2/22/2008
KB944533	Cumulative Security Update for Internet Explorer 6 Service Pack 1	Yes	Pass	2/22/2008
KB890830	Windows Malicious Software Removal Tool - January 2008	Yes	Pass	1/14/2008
KB941644	Security Update for Windows 2000	Yes	Pass	1/14/2008
KB943485	Security Update for Windows 2000	Yes	Pass	1/14/2008
KB958687	Security Update for Windows 2000	Yes	Pass	1/14/2008
KB890830	Windows Malicious Software Removal Tool - December 2007	Yes	Pass	12/17/2007
KB937894	Security Update for Windows 2000	Yes	Pass	12/17/2007
KB941568	Security Update for DirectX 9.0 for Windows 2000	Yes	Pass	12/17/2007
KB941569	Security Update for Windows 2000	Yes	Pass	12/17/2007
KB942615	Cumulative Security Update for Internet Explorer 6 Service Pack 1	Yes	Pass	12/17/2007
KB890830	Windows Malicious Software Removal Tool - November 2007	Yes	Pass	11/29/2007
KB890830	Windows Malicious Software Removal Tool - October 2007	Yes	Pass	10/15/2007
KB923810	Security Update for Windows 2000	Yes	Pass	10/15/2007
KB933729	Security Update for Windows 2000	Yes	Pass	10/15/2007
KB939653	Cumulative Security Update for Internet Explorer 6 Service Pack 1	Yes	Pass	10/15/2007
KB941202	Security Update for Outlook Express 6 Service Pack 1	Yes	Pass	10/15/2007
KB890830	Windows Malicious Software Removal Tool - September 2007	Yes	Pass	9/22/2007
KB938827	Security Update for Windows 2000	Yes	Pass	9/22/2007
KB890830	Windows Malicious Software Removal Tool - August 2007	Yes	Pass	8/22/2007
KB921503	Security Update for Windows 2000	Yes	Pass	8/22/2007
KB936021	Security Update for Windows 2000	Yes	Pass	8/22/2007
KB936782	Security Update for Windows Media Player 9	Yes	Pass	8/22/2007
KB937143	Cumulative Security Update for Internet Explorer 6 Service Pack 1	Yes	Pass	8/22/2007
KB938127	Security Update for Internet Explorer 6 Service Pack 1	Yes	Pass	8/22/2007
KB938829	Security Update for Windows 2000	Yes	Pass	8/22/2007
KB890830	Windows Malicious Software Removal Tool - July 2007	Yes	Pass	7/16/2007
KB926122	Security Update for Windows 2000	Yes	Pass	7/16/2007
KB928366	Security Update for Microsoft .NET Framework, Version 1.1, Service Pack 1	Yes	Pass	7/16/2007
KB890830	Windows Malicious Software Removal Tool - June 2007	Yes	Pass	6/19/2007
KB927891	Update for Windows 2000	Yes	Pass	6/19/2007
KB933566	Cumulative Security Update for Internet Explorer 6 Service Pack 1	Yes	Pass	6/19/2007
KB935839	Security Update for Windows 2000	Yes	Pass	6/19/2007
KB935840	Security Update for Windows 2000	Yes	Pass	6/19/2007
KB890830	Windows Malicious Software Removal Tool - May 2007	Yes	Pass	5/21/2007
KB890830	Windows Malicious Software Removal Tool - April, 2007	Yes	Pass	5/21/2007
KB931768	Cumulative Security Update for Internet Explorer 6 Service Pack 1	Yes	Pass	5/21/2007
KB931906	Security Update for CAPICOM	Yes	Pass	5/21/2007
KB925902	Security Update for Windows 2000	Yes	Pass	4/13/2007
KB930178	Security Update for Windows 2000	Yes	Pass	4/13/2007
KB931784	Security Update for Windows 2000	Yes	Pass	4/13/2007
KB932168	Security Update for Windows 2000	Yes	Pass	4/13/2007
KB890830	Windows Malicious Software Removal Tool - March 2007	Yes	Pass	3/21/2007
KB890830	Windows Malicious Software Removal Tool - February 2007	Yes	Pass	2/16/2007
KB918118	Security Update for Windows 2000	Yes	Pass	2/16/2007
KB924667	Security Update for Window 2000	Yes	Pass	2/16/2007

KB926436	Security Update for Window 2000	Yes	Pass	2/16/2007
KB927779	Security Update for Microsoft Data Access Components 2.5 Service Pack 3	Yes	Pass	2/16/2007
KB928090	Cumulative Security Update for Internet Explorer 6 Service Pack 1	Yes	Pass	2/16/2007
KB928843	Security Update for Window 2000	Yes	Pass	2/16/2007
KB890830	Windows Malicious Software Removal Tool - January 2007	Yes	Pass	1/15/2007
KB890830	Windows Malicious Software Removal Tool - Dec 2006	Yes	Pass	12/14/2006
KB923689	Security Update for Window 2000	Yes	Pass	12/14/2006
KB923694	Cumulative Security Update for Outlook Express 6 SP1	Yes	Pass	12/14/2006
KB925398	Security Update for Windows Media Player 6.4	Yes	Pass	12/14/2006
KB925454	Cumulative Security Update for IE 6 SP 1 (MS Bulletin MS06-072)	Yes	Pass	12/14/2006
KB890830	Windows Malicious Software Removal Tool - November 2006	Yes	Pass	11/20/2006
KB920213	Security Update for Windows 2000	Yes	Pass	11/20/2006
KB922760	Cumulative Security Update for Internet Explorer 6 Service Pack1	Yes	Pass	11/20/2006
KB923980	Security Update for Windows 2000	Yes	Pass	11/20/2006
KB924270	Security Update for Windows 2000	Yes	Pass	11/20/2006
KB890830	Windows Malicious Software Removal Tool - October 2006	Yes	Pass	10/19/2006
KB920958	Security Update for Windows 2000	Yes	Pass	10/19/2006
KB923191	Security Update for Windows 2000	Yes	Pass	10/19/2006
KB923272	Security Update for Office 2003	Yes	Pass	10/19/2006
KB923414	Security Update for Windows 2000	Yes	Pass	10/19/2006
KB924191	Security Update for Windows 2000	Yes	Pass	10/19/2006
KB925486	Security for Internet Explorer 6 Service Pack 1	Yes	Pass	10/19/2006
KB20683	Security Update for Windows 2000	Yes	Pass	9/15/2006
KB917008	Security Update for Windows 2000	Yes	Pass	9/15/2006
KB917422	Security Update for Windows 2000	Yes	Pass	9/15/2006
KB918899	Cumulative Security Update for IE 6 SP1	Yes	Pass	9/15/2006
KB920670	Security Update for Windows 2000	Yes	Pass	9/15/2006
KB920958	Security Update for Windows 2000	Yes	Pass	9/15/2006
KB921398	Security Update for Windows 2000	Yes	Pass	9/15/2006
KB921883	Security Update for Windows 2000	Yes	Pass	9/15/2006
KB922582	Update for Windows 2000	Yes	Pass	9/15/2006
KB922612	Security Update for Windows 2000	Yes	Pass	9/15/2006
KB973540	Security Update for Windows Media Player 9 for windows 2000	Yes	Pass	8/12/2006
KB917537	Security Update for Windows 2000	Yes	Pass	8/3/2006
KB914388	Security Update for Windows 2000	Yes	Pass	7/14/2006
KB917159	Security Update for Windows 2000	Yes	Pass	7/14/2006
KB911280	Security Update for Windows 2000	Yes	Pass	6/13/2006
KB913580	Security Update for Windows 2000	Yes	Pass	6/13/2006
KB914389	Security Update for Windows 2000	Yes	Pass	6/13/2006
KB916281	Cumulative Security Update for Internet Explorer 6 Service Pack 1	Yes	Pass	6/13/2006
KB917344	Security Update for (Jscript Version 5.6) for windows 2000	Yes	Pass	6/13/2006
KB917734	Security Update for Windows Media Player 9	Yes	Pass	6/13/2006
KB917736	Security Update for Windows 2000	Yes	Pass	6/13/2006
KB917953	Security Update for Windows 2000	Yes	Pass	6/13/2006
KB908531	Security Update for Windows 2000	Yes	Pass	4/18/2006
KB911562	Security Update for IE 6 SP1 Windows 2000	Yes	Pass	4/18/2006
KB911564	Security Update for Windows Media Player Plug-in	Yes	Pass	4/18/2006
KB911565	Security Update for Windows Media Player 9	Yes	Pass	4/18/2006
KB911567	Cumulative Security Update for Outlook Express 6 SP1	Yes	Pass	4/18/2006
MSXML 4.0	Security Update for Windows 2000	Yes	Pass	4/18/2006
KB890830	Windows Malicious Software Removal Tool - January 2006	Yes	Pass	1/10/2006
KB892130	Windows Genuine Advantage Validation Tool	Yes	Pass	1/10/2006
KB896423	Security Update for Windows 2000	Yes	Pass	1/10/2006
KB896424	Security Update for Windows 2000	Yes	Pass	1/10/2006
KB899589	Security Update for Windows 2000	Yes	Pass	1/10/2006
KB900725	Security Update for Windows 2000	Yes	Pass	1/10/2006
KB901017	Security Update for Windows 2000	Yes	Pass	1/10/2006
KB902400	Security Update for Windows 2000	Yes	Pass	1/10/2006
KB904706	Update for DirectX 9 for Windows 2000	Yes	Pass	1/10/2006
KB905414	Security Update for Windows 2000	Yes	Pass	1/10/2006
KB905495	Security Update for IE 6 SP1 Windows 2000	Yes	Pass	1/10/2006
KB905749	Security Update for Windows 2000	Yes	Pass	1/10/2006
KB905915	Cumulative Security Update for IE 6 SP1	Yes	Pass	1/10/2006
KB908519	Security Update for Windows 2000	Yes	Pass	1/10/2006

KB908523	Security Update for Windows 2000	Yes	Pass	1/10/2006
KB912919	Security Update for Windows 2000	Yes	Pass	1/10/2006
* Also Tested	Microsoft AntiSpyware Beta 1	Yes	Pass	8/29/2005
KB890046	Security Update for Windows 2000	Yes	Pass	8/29/2005
KB891861	Update Rollup 1 for Windows 2000	Yes	Pass	8/29/2005
KB893756	Security Update for Windows 2000	Yes	Pass	8/29/2005
KB896243	Security Update for Windows 2000 - Addresses a critical Security problem	Yes	Pass	8/29/2005
KB896358	Security Update for Windows 2000	Yes	Pass	8/29/2005
KB896422	Security Update for Windows 2000	Yes	Pass	8/29/2005
KB896727	Cumulative Security Update for Internet Explorer 6 SP1	Yes	Pass	8/29/2005
KB897715	Security Update for Outlook Express 6 SP 1	Yes	Pass	8/29/2005
KB899587	Security Update for Windows 2000	Yes	Pass	8/29/2005
KB899588	Security Update for Windows 2000 - Addresses a critical Security problem	Yes	Pass	8/29/2005
KB901214	Security Update for Windows 2000	Yes	Pass	8/29/2005
KB893086	Security Update for Windows 2000	Yes	Pass	5/24/2005
KB894320	Security Update for Windows 2000	Yes	Pass	5/24/2005
KB831167	Critical Update for Internet Explorer 6 Service Pack 1	Yes	Pass	5/21/2005
KB842773	BITS 2.0 and Win HTTP 5.1	Yes	Pass	4/27/2005
KB890859	Security Update for Windows 2000	Yes	Pass	4/27/2005
KB890923	Cumulative Security Update for Internet Explorer 6 Service Pack 1	Yes	Pass	4/27/2005
KB893066	Security Update for Windows 2000	Yes	Pass	4/27/2005
KB893096	Security Update for Windows 2000	Yes	Pass	4/27/2005
n/a	Microsoft Windows Installer 3.1	Yes	Pass	4/27/2005
KB867460	Microsoft .Net Framework 1.1 Service Pack 1	Yes	Pass	3/8/2005
KB886903	Security Update, Microsoft .Net Framework v1.1 Service Pack 1	Yes	Pass	3/8/2005
330994	April 2003, Security Update for Outlook Express 6 SP1	Yes	Pass	2/14/2005
814078	Security Update (Microsoft Jscript version 5.6, Windows 2000, Windows XP)	Yes	Pass	2/14/2005
814078	Security Update (Microsoft Jscript version 5.1, Windows 2000) Web site	Yes	Pass	2/14/2005
818043	Recommended Update for Windows 2000 Web site	Yes	Pass	2/14/2005
822831	Recommended Update for Windows 2000	Yes	Pass	2/14/2005
822831	Recommended Update for Windows 2000	Yes	Pass	2/14/2005
823559	Security Update for Microsoft Windows	Yes	Pass	2/14/2005
823559	Security Update for Microsoft Windows Web site	Yes	Pass	2/14/2005
KB329115	Windows 2000 Hotfix	Yes	Pass	2/14/2005
KB820888	Update for Windows 2000	Yes	Pass	2/14/2005
KB822831	Windows 2000 Hotfix	Yes	Pass	2/14/2005
KB823182	Security Update for Windows 2000	Yes	Pass	2/14/2005
KB823353	Cumulative Security Update for Outlook Express 6 SP1	Yes	Pass	2/14/2005
KB823559	Windows 2000 Hotfix	Yes	Pass	2/14/2005
KB823980	Windows 2000 Hotfix	Yes	Pass	2/14/2005
KB824105	Security Update for Microsoft Windows	Yes	Pass	2/14/2005
KB824141	Security Update for Microsoft Windows	Yes	Pass	2/14/2005
KB824141	Windows 2000 Hotfix	Yes	Pass	2/14/2005
KB824146	Security Update for Windows 2000	Yes	Pass	2/14/2005
KB824146	Windows 2000 Hotfix	Yes	Pass	2/14/2005
KB825119	Security Update for Microsoft Windows 2000	Yes	Pass	2/14/2005
KB825119	Windows 2000 Hotfix	Yes	Pass	2/14/2005
KB826232	Security Update for Microsoft Windows 2000	Yes	Pass	2/14/2005
KB826232	Windows 2000 Hotfix	Yes	Pass	2/14/2005
KB828026	Critical Update for Windows Media Player Script Commands	Yes	Pass	2/14/2005
KB828028	Security Update for Windows 2000	Yes	Pass	2/14/2005
KB828028	Windows 2000 Hotfix	Yes	Pass	2/14/2005
KB828028	Windows 2000 Hotfix	Yes	Pass	2/14/2005
KB828035	Security Update for Microsoft Windows 2000	Yes	Pass	2/14/2005
KB828035	Windows 2000 Hotfix	Yes	Pass	2/14/2005
KB828035	Windows 2000 Hotfix	Yes	Pass	2/14/2005
KB828741	Security Update for Windows 2000	Yes	Pass	2/14/2005
KB828741	Windows 2000 Hotfix	Yes	Pass	2/14/2005
KB828749	Security Update for Microsoft Windows	Yes	Pass	2/14/2005
KB828749	Windows 2000 Hotfix	Yes	Pass	2/14/2005
KB829558	Update for Jet 4.0 Service Pack 8	Yes	Pass	2/14/2005
KB832483	Security Update for Microsoft Data Access Components	Yes	Pass	2/14/2005
KB832894	Cumulative Security Update for Internet Explorer 6 Service Pack 1	Yes	Pass	2/14/2005
KB832894	Cumulative Security Update for Internet Explorer 5.01 Service Pack 4	Yes	Pass	2/14/2005

KB833989	Security Update for Internet Explorer 6 Service Pack 1	Yes	Pass	2/14/2005
KB834707	Cumulative Security Update for Internet Explorer 6 Service Pack 1	Yes	Pass	2/14/2005
KB835732	Security Update for Windows 2000	Yes	Pass	2/14/2005
KB837001	Security Update for Windows 2000	Yes	Pass	2/14/2005
KB837272	Update for Windows Media Player 9 Series	Yes	Pass	2/14/2005
KB839643	Security Update for DirectX 9.0	Yes	Pass	2/14/2005
KB839645	Security Update for Windows 2000	Yes	Pass	2/14/2005
KB839645	Windows 2000 Hotfix	Yes	Pass	2/14/2005
KB840315	Security Update for Windows 2000	Yes	Pass	2/14/2005
KB840315	Windows 2000 Hotfix	Yes	Pass	2/14/2005
KB840987	Security Update for Windows 2000	Yes	Pass	2/14/2005
KB840987	Windows 2000 Hotfix	Yes	Pass	2/14/2005
KB841356	Security Update for Windows 2000	Yes	Pass	2/14/2005
KB841356	Windows 2000 Hotfix	Yes	Pass	2/14/2005
KB841533	Security Update for Windows 2000	Yes	Pass	2/14/2005
KB841533	Windows 2000 Hotfix	Yes	Pass	2/14/2005
KB841872	Security Update for Windows 2000	Yes	Pass	2/14/2005
KB841872	Windows 2000 Hotfix	Yes	Pass	2/14/2005
KB841873	Security Update for Windows 2000	Yes	Pass	2/14/2005
KB841873	Windows 2000 Hotfix	Yes	Pass	2/14/2005
KB842526	Security Update for Windows 2000	Yes	Pass	2/14/2005
KB867282	Security Update - Internet Explorer 6	Yes	Pass	2/14/2005
KB867801	Cumulative Security Update for Internet Explorer 6 Service Pack 1	Yes	Pass	2/14/2005
KB870669	Critical Update for ADODB.stream	Yes	Pass	2/14/2005
KB871250	Security Update for Windows 2000	Yes	Pass	2/14/2005
KB873333	Security Update for Windows 2000	Yes	Pass	2/14/2005
KB873339	Security Update for Windows 2000	Yes	Pass	2/14/2005
KB873374	Microsoft GDI+ Detection Tool	Yes	Pass	2/14/2005
KB885250	Security Update for Windows 2000	Yes	Pass	2/14/2005
KB885492	Security Update for Windows Media Player 9 Series	Yes	Pass	2/14/2005
KB885835	Security Update for Windows 2000	Yes	Pass	2/14/2005
KB885835	Windows 2000 Hotfix	Yes	Pass	2/14/2005
KB885836	Security Update for Windows 2000	Yes	Pass	2/14/2005
KB885836	Windows 2000 Hotfix	Yes	Pass	2/14/2005
KB886179	Security Update for Windows Journal Viewer	Yes	Pass	2/14/2005
KB887797	Cumulative Update for Outlook Express 6 Service Pack 1	Yes	Pass	2/14/2005
KB888113	Security Update for Windows 2000	Yes	Pass	2/14/2005
KB889293	Cumulative Security Update for Internet Explorer 6 Service Pack 1	Yes	Pass	2/14/2005
KB890047	Security Update for Windows 2000	Yes	Pass	2/14/2005
KB890175	Security Update for Windows 2000	Yes	Pass	2/14/2005
KB891711	Security Update for Windows 2000	Yes	Pass	2/14/2005
KB891781	Security Update for Windows 2000	Yes	Pass	2/14/2005
n/a	Microsoft Internet Explorer 6 Service Pack 1 (Windows 2000)	Yes	Pass	2/14/2005
n/a	Microsoft .NET Framework Version 1.1	Yes	Pass	2/14/2005
Q329115	Security Update (Windows 2000)	Yes	Pass	2/14/2005
Q818043	Windows 2000 Hotfix (SP5)	Yes	Pass	2/14/2005
WM828026	Windows Media Player Hotfix	Yes	Pass	2/14/2005
KB929969	Security Update for Internet Explorer 6 Service Pack 1	Yes	Pass	1/15/2005