



Your Vision, Our Future

March 14, 2012

The following list of Lenovo provided system updates has been tested for **EndoWorks 7.4** compatibility. Prior to applying System Driver or BIOS updates, make sure your system is current with EndoWorks patches. Olympus cannot guarantee the current content of the OEM website; consequently, the specific approved update may no longer be available. If this is the case, please do not install a different update which has not yet been tested. For questions regarding applying System updates or EndoWorks patches, please contact Olympus' Technical Assistance Center (TAC) at 1-800-848-9024.

Lenovo system updates are available at http://support.lenovo.com/en_US/downloads/default.page

In Type select Desktops & all-in-ones, in Series select ThinkCentre M Series desktops, in Subseries select ThinkCentre M57 or M58p, then in Machine Type select as listed below

| M57 Machine Type/Model | BIOS/Driver Description | Operating System | Version and Release Date | Tested | Pass/Fail | Date |
|--------------------------------|--|--|-------------------------------------|---------------|------------------|-------------|
| 6075-BJU | BIOS | Microsoft Windows XP Professional (32 Bit) | 62, 7/19/2011 | Yes | Pass | 3/14/2012 |
| 6075-BJU | Intel 82566DM-2 Gigabit NIC | Microsoft Windows XP Professional (32 Bit) | 9.13.16.0, 2/2/2010 | Yes | Pass | 8/25/2011 |
| 6075-BJU | Intel Q35 express chipset Display Adapter | Microsoft Windows XP Professional (32 Bit) | 6.14.10.5218, 1/13/2010 | Yes | Pass | 8/25/2011 |
| 6075-BJU | ADI Sound MAX Audio | Microsoft Windows XP Professional (32 Bit) | 5.10.01.6380, 6/22/2007 | Yes | Pass | 8/25/2011 |
| M58p Machine Type/Model | BIOS/Driver Description | Operating System | Description and Release Date | Tested | Pass/Fail | Date |
| 6138-A1U or 6138-G71 | BIOS | Microsoft Windows XP Professional (32 Bit) | 71, 7/27/2011 | Yes | Pass | 3/14/2012 |
| 6138-A1U or 6138-G71 | Intel® 8256LM-3 Gigabit NIC | Microsoft Windows XP Professional (32 Bit) | 11.2.19.0, 9/23/2009 | Yes | Pass | 8/25/2011 |
| 6138-A1U or 6138-G71 | Intel® Q45/Q43 Express Chipset Display Adapter | Microsoft Windows XP Professional (32 Bit) | 6.14.10.5212, 1/8/2010 | Yes | Pass | 8/25/2011 |
| 6138-A1U or 6138-G71 | ADI Sound MAX Audio | Microsoft Windows XP Professional (32 Bit) | 5.10.01.6540-0005, 4/30/2009 | Yes | Pass | 8/25/2011 |



Your Vision, Our Future

March 14, 2012

The following list of IBM provided system updates has been tested for **EndoWorks 7.4** compatibility. Prior to applying System Driver or BIOS / UEFI updates, make sure your system is current with EndoWorks patches. Olympus cannot guarantee the current content of the OEM website; consequently, the specific approved update may no longer be available. If this is the case, please do not install a different update which has not yet been tested. For questions regarding applying System updates or EndoWorks patches, please contact Olympus' Technical Assistance Center (TAC) at 1-800-848-9024.

IBM system updates are available at <http://www-947.ibm.com/support/entry/portal/Downloads>

Under "1. Choose your products" select the "Search for a product" radio button, then scroll down and click Continue. Under "Downloads and fixes", select the Search for downloads and fixes" link. Type either x3400 or x3400 M3 in the Search for text box, then click the Search icon next to the field. In the displayed search results, select the appropriate System Driver or BIOS /UEFI update link, then wait while the webpage is refreshed with the appropriate file links.

| x3400 Product | BIOS/Driver Description | Operating System | Version and Release Date | Tested | Pass/Fail | Date |
|-------------------------|-------------------------------------|-------------------------------|-------------------------------------|---------------|------------------|-------------|
| 7976 | BIOS | Microsoft Windows Server 2003 | 1.74, 11/25/2010 | Yes | Pass | 8/25/2011 |
| 7976 | RAID Firmware | Microsoft Windows Server 2003 | 5.2.0-17002, 07/27/2011 | Yes | Pass | 3/14/2012 |
| 7976 | Broadcom NetXtreme Gigabit Ethernet | Microsoft Windows Server 2003 | 10.78.0.0, 3/19/2008 | Yes | Pass | 8/25/2011 |
| 7976 | ATI ES1000 Display Adapter | Microsoft Windows Server 2003 | 8.24.50.2000, 10/25/2007 | Yes | Pass | 8/25/2011 |
| x3400 M3 Product | BIOS/Driver Description | Operating System | Description and Release Date | Tested | Pass/Fail | Date |
| 7379 | UEFI | Microsoft Windows Server 2003 | 1.09, 6/29/2011 | Yes | Pass | 3/14/2012 |
| 7379 | RAID Firmware | Microsoft Windows Server 2003 | 11.0.1-0042, 7/28/2011 | Yes | Pass | 3/14/2012 |
| 7379 | Broadcom NetXtreme Gigabit Ethernet | Microsoft Windows Server 2003 | 6.0.32.0, 11/9/2010 | Yes | Pass | 3/14/2012 |