

Biopsy Forceps

FB-56D-1

Disposable Cytology Brush

BC-201C-1006/BC-203D-2006

Disposable Grasping Forceps

FG-51D/52D/54D/55D



OLYMPUS ENDO-THERAPY

Ideal for Diagnosis of Peripheral Regions



Now you can take advantage of bronchoscopic instrumentation in previously difficult to approach peripheral regions.

Biopsy Forceps



FB-56D-1

FB-56D-1

Featuring a rat tooth jaw, these reusable forceps provide greater precision and reliability for tissue collection. The unique Olympus two stage flexibility at the distal end offers improved insertion capability when the scope tip is angulated.

Working length : 1,150mm

Disposable Cytology Brush



BC-201C-1006



BC-203D-2006

BC-201C-1006

Contains a kink-resistant wire to improve insertion capability. For your convenience, the single use brush is packed sterile, ready to use.

Working length : 1,050mm
Brush diameter : 1mm
Brush length : 6mm


BC-203D-2006

Features a covered sheath that improves specimen collection and prevents cross-contamination. For your convenience, the single use brush is packed sterile, ready to use.

Working length : 1,150mm
Brush diameter : 2mm
Brush length : 6mm



FG-51D



FG-52D



FG-54D



FG-55D

FG-51D/52D/54D/55D

Contains a smooth sheath that facilitates passage through the scope channel. The baskets have been specially engineered allowing them to be opened widely, making it well suited for retrieval of foreign objects. For your convenience, the single use forceps is packed sterile, ready to use.

FG-51D
Working length : 1,150mm
Spiral basket type with 4 spiral wires
Basket width : 11mm

FG-52D
Working length : 1,150mm
Spiral basket type with 3 spiral wires
Basket width : 9mm

FG-54D
Working length : 1,150mm
Three-claw type
Claw width : 10mm

FG-55D
Working length : 1,150mm
Basket type with 4 parallel wires
Basket width : 14mm



Handle Section

Specification, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

OLYMPUS

OLYMPUS OPTICAL CO., LTD.
2951 Ishikawa-cho, Hachioji-shi, Tokyo 192-8507, Japan
OLYMPUS OPTICAL CO. (EUROPA) GMBH.
Postfach 10 49 08, 20034 Hamburg / Wendenstrasse 14-18, 20097 Hamburg, Germany
OLYMPUS AMERICA INC.
2 Corporate Center Drive, Melville, N.Y. 11747-3157, U.S.A.
OLYMPUS LATIN AMERICA, INC.
6100 Blue Lagoon Drive, Suite 390 Miami, Florida 33126-2087, U.S.A
OLYMPUS KEYMED
Keymed House, Stock Road, Southend-on-Sea, Essex SS2 5QH, England

OLYMPUS SINGAPORE PTE LTD.
491B River Valley Road #12-01/04, Valley Point Office Tower, Singapore 248373
OLYMPUS HONG KONG AND CHINA LIMITED.
Room 1520-1527, Ocean Centre, 5 Canton Road, Tsimshatsui, Kowloon, Hong Kong
OLYMPUS BEIJING REPRESENTATIVE OFFICE
Room 401, East Tower, Gong Yuan No6 Royal Palace, No6 Gong YuanXijie, JianGuoMenNei, DongCheng District, Beijing, 100005, China
OLYMPUS MOSCOW LIMITED LIABILITY COMPANY
117071, Moscow, Malaya Kaluzhskaya 19, bld. 1, fl.2, Russia
OLYMPUS AUSTRALIA PTY, LTD.
31 Gilby Road, Mount Waverley, VIC., 3149, Australia

www.olympus.com

Printed in Japan F1002SB-0303