

# OLYMPUS

Your Vision, Our Future

## Olympus Endoscopic Mouthpiece

Single use endoscopic mouthpiece with rear lip guard MAJ-168



**Assured scope compatibility and handling.  
Optimized patient comfort.**

OLYMPUS

Endo**Therapy**

# OLYMPUS QUALITY, DESIGN AND COMPATIBILITY

The Olympus Endoscopic Mouthpiece meets today's needs of a busy GI suite. An increasing number of procedures are being performed and the pressure for fast diagnosis and treatment of patients is a major factor in creating a successful endoscopy department.

Olympus MAJ-168 Single Use Mouthpiece was designed to optimize patient comfort and care of gums and teeth. Ergonomically shaped rear lip guard prevents trauma, injury and scope damage. Provided non-latex strap prevents accidental mouthpiece shifting.

Olympus MAJ-168 Single Use Mouthpiece is perfectly fitted and sized to fit all your Olympus scopes.

- Side connectors and an adjustable strap allow the mouthpiece to be fitted comfortably to the patient.
- Silicone rubber retaining strap prevents the risk of allergy from latex based equivalents.
- Supplied clean in individually wrapped in boxes of 50 mouthpieces.
- Packaging size kept to a minimum with access either through a front flap or from the top of the box, depending on user preference.

## SPECIFICATION

The Olympus mouthpiece is manufactured in the UK from polyethylene and is provided with a silicone rubber retaining strap, both manufactured to meet the requirements of ISO 10993 and USP CLASS VI. Aperture is obround in shape and is a minimum of 20mm (H) x 27mm (W).

Olympus reserves the right to alter the specification without notice.

# OLYMPUS®

**OLYMPUS CORPORATION**  
Shinjuku Monolith, 3-1, Nishi Shinjuku 2-chome,  
Shinjuku-ku, Tokyo 163-0914, Japan  
**OLYMPUS KEYMED**  
KeyMed House, Stock Road, Southend-on-Sea,  
Essex SS2 5QH, United Kingdom

**OLYMPUS EUROPA GMBH**  
(Premises/Goods delivery) Wendenstraße 14-18, 20097 Hamburg, Germany  
(Letters) Postfach 10 49 08, D-20034 Hamburg, Germany  
**OLYMPUS AMERICA INC.**  
Two Corporate Center Drive, Melville,  
New York 11747-3157, USA  
**OLYMPUS SINGAPORE PTE LTD.**  
491B River Valley Road #12-01/04, Valley Point Office Tower,  
Singapore 248373, Singapore

KeyMed Limited



ISO 9001 Registered  
Certificate No. FM 20993  
Quality Management System



ISO 14001 Registered  
Certificate No. EMS 65964  
Environmental Management  
System