

OLYMPUS®



EVIS EXERA GASTROINTESTINAL VIDEOSCOPE
OLYMPUS GIF TYPE Q160


LARGE DISPLAY WITH SUPERIOR IMAGE QUALITY

Large images with superior sharpness and clarity, expanded field of view, and great depth of field make this high-quality scope the right choice for a wide variety of gastroscopic examination and treatment procedures





Main Features for Superb Performance

- Large display size delivers superior image quality with sharper, crisper image and true-to-life color from edge to edge.
- Extra-wide 140° field of view enables accurate observation of a wider area.
- Slim 9.5 mm outer diameter insertion tube with 2.8 mm channel diameter.
- 4-way angulation (210° up, 90° down and 100° right/left) enables complete and comprehensive examination of the upper digestive tract.
- Ergonomically designed grip enhances scope maneuverability while redesigned easy-to-access knobs and user-selectable switches improve operability.
- Fully compatible with the CV-160/140.
- Scope ID function stores individual scope information and displays it on the monitor to facilitate endoscopy suite management.





Specifications

Optical System	Field of view	140 _i
	Direction of view	0 _j Forward viewing
	Depth of field	3 to 100 mm
Distal End	Outer diameter	9.8 mm
	Insertion Tube	Outer diameter
Bending Section	Angulation range	Up 210 _i , Down 90 _i , Right 100 _i , Left 100 _i
Working Length		1030 mm
Total Length		1345 mm
Instrument Channel	Inner diameter	2.8 mm
	Minimum visible distance	3 mm from the distal end
	Endo-Therapy accessory entrance/exit position in field of view 	



Olympus business areas



OLYMPUS
FOCUS ON LIFE

<http://www.olympus.com>

OLYMPUS OPTICAL CO., LTD.

San-Ei Building, 22-2, Nishi Shinjuku 1-chome, Shinjuku-ku, Tokyo, Japan

OLYMPUS OPTICAL CO. (EUROPA) GMBH,
Postfach 10 49 08, 20034 Hamburg / Wendenstrasse 14-16, 20097 Hamburg, Germany

OLYMPUS AMERICA INC.,
2 Corporate Center Drive, Melville, N.Y. 11747-3157, U.S.A.

OLYMPUS KEYMED
Keymed House, Stock Road, Southend-on-Sea, Essex SS2 5QH, England

OLYMPUS SINGAPORE PTE LTD.
491B River Valley Road #12-01/04 Valley Point Office Tower, Singapore 248373

OLYMPUS BEIJING REPRESENTATIVE OFFICE
Rm. 818, South Tower, Beijing Kerry Centre, No.1 Guanghua Road, Chaoyang District, Beijing, 100020, China

OLYMPUS MOSCOW LIMITED LIABILITY COMPANY
117071, Moscow, Malaya Kaluzhskaya 19, bld. 1, fl.2, Russia

OLYMPUS AUSTRALIA PTY. LTD.
1/104 Ferntree Gully Road, Oakleigh, VIC 3166, Australia