TC-EU COMPACT ENDOSCOPIC ULTRASOUND WORKSTATION

A DEDICATED COMPACT ENDOSCOPIC ULTRASOUND WORKSTATION

The Olympus Endoscopy Support range of products includes a complete family of workstations, developed to be practical and space-efficient wherever they are used. Ergonomically designed, with many innovative features, they are more than just trolleys, they are unique, purpose-designed mobile workstations.

The Olympus TC-EU is a compact workstation designed to support the EU-M30S Endoscopic Ultrasound processor.

FEATURES INCLUDE:

- **Compact** space-efficient for use in what is so often a cramped working environment.
- Aesthetically pleasing design to suit the modern clinical environment.
- Base tray designed to accommodate ancillary items for the EU-M30S.
- Mobile with four 100mm castors, two diagonally opposed with brakes. Easily brought into or moved out of the procedure area.
- Colour co-ordinated bumper protection on all four corners - to prevent accidental damage to walls, doors and paintwork.
- Keyboard shelf easy keyboard access for operation of EU-M30S processor.
- **Top shelf** to accommodate a 14" monitor.
- **Fixed second shelf** to accommodate the EU-M30S processor.
- Adjustable third shelf to accommodate user specified equipment, such as video printer or video recorder.
- Integral handle ideally positioned for easy manoeuvrability.
- **Back Panel** is available to provide a fully enclosed workstation.
- Plug Blocks a multiway plug block kit is available reducing the number of trailing cables

CONSTRUCTION:

Olympus workstations, designed and manufactured by KeyMed, have an enviable reputation for performance and durability the TC-EU is no exception.

Built to an exhaustive quality standard, the all metal frame and shelf construction are electrostatically powder painted and will be highly durable within the operating environment.





SPECIFICATION:

Dimensions: $480(W) \times 1150(H) \times 570(D) \text{ mm}$

Unladen Weight: 29 kg

Maximum

carrying capacity: 100 kg

Olympus reserves the right to alter the specification without notice.



Issue 1: September 2002

OLYMPUS OPTICAL CO., LTD. San-Ei Building, 22-2, Nishi Shinjuku 1-chome, Shinjuku-ku, Tokyo, Japan OLYMPUS KEYMED OLYMPUS OPTICAL CO., (EUROPA) GMBH.

(Premises/Goods delivery) Wendenstra§e 14-18, 20097 Hamburg, German

OLYMPUS AMERICA INC.

wo Corporate Center Drive, Melville, New York 11747-3157, USA

OLYMPUS SINGAPORE PTE LTD.

Printed in UK KI 5047994/902



OLYMPUS

TC-EU COMPACT ENDOSCOPIC ULTRASOUND WORKSTATION



A further product in a design co-ordinated endoscopy support range

