

PRODUCT COMPARISON CHART

OLYMPUS[®]

Your Vision, Our Future

EBUS Bronchoscopes	Universal Platform	NBI	HD	Wide Angle (170)	Close Focus (2mm)
BF-UC160F-OL8	x				
Fiber-optic Bronchoscopes					
BF-XP60	•				
BF-MP60	•				
BF-P60	•				
BF-1T60	•				
BF-3C40	•				
BF-XT40	•				
Video Bronchoscopes					
BF-Q180	x	x			
BF-Q180-AC	x	x			
BF-P180	x	x			
BF-1T180	x	x			
BF-1TQ180	x	x			
LTF-160	x				
BF-XP160F	x				
BF-MP160F	x				
BF-3C160	x				
BF-XT160	x				
Capsule Endoscopes					
Endo Capsule	x				
Choledochoscopes					
CHF-BP30	•				
Colonoscopes:					
CF-H180AL/I	x	x	x	x	x
CF-Q180AL/I	x	x		x	
PCF-Q180AL/I	x	x			
CF-Q160ZL/I	x	x			
PCF-H180AL/I	x	x	x		x
CF-2T160L/I	x				
Duodenoscopes					
TJF-Q180V	x	x			
Echoendoscopes					
GF-UE160-AL5	x				
GF-UM160	x				
GF-UC140P-AL5 / GF-UCT140-AL5	x				
GF-UC160P-OL5 / GF-UCT160-OL5	x				

• **Fiberscopes - can connect to the CV-180 with coupler adaptor and camera head**

PRODUCT COMPARISON CHART



Your Vision, Our Future

	Universal Platform	NBI	HD	Wide Angle (170)	Close Focus (2mm)
Enteroscopes					
SIF-Q180	x	x			
Esophagoscopes					
PEF-V	x				
Gastrosopes					
GIF-H180J	x	x	x		x
GIF-H180	x	x	x		x
GIF-Q180	x	x			
GIF-N180	x	x			
GIF-XP180N	x	x			
GIF-1TQ160	x				
GIF-2T160	x				
GIF-XTQ160	x				
GIF-Q160Z	x	x			
Laryngoscopes					
LF-GP	•				
LF-TP	•				
LF-DP	•				
Rhinolaryngoscopes					
ENF-V2	x	x			
ENF-VQ	x	x			
ENF-VT2	x				x
ENF-GP	•				
ENF-P4	•				
ENF-XP	•				
Sigmoidoscopes					
PCF-S	x				
CF-Q160S	x				
OSF-V60					
OSF-4					
• Fiberscopes - can connect to the CV-180 with coupler adaptor and camera head					